

FOREHEAD THERMOMETER CONTACT FREE 4

apornorm[®]
die marke der apotheke

Non-contact measurement



Distance
control



Non-contact
measurement

1 sec

Measurement
time





**CLINICALLY +
PROVEN**

5 year
warranty



A quality
brand of

WEPA
DIE APOTHEKENMARKE



Switches & housings

- ① Measuring probe
- ② Measuring button
- ③ SET button
- ④ Battery compartment cover
- ⑤ Display

Display

- ⑥ Object mode
- ⑦ Body mode
- ⑧ Distance indicator
- ⑨ Sound on/off
- ⑩ Battery level
- ⑪ Celsius (°C) or Fahrenheit (°F) unit display
- ⑫ Memory mode
- ⑬ Hour mode (AM/PM)
- ⑭ Date/time display
- ⑮ Temperature value

Buttons	Functions
Measuring button	Start measurement/switch off device
SET button	Sound On/Off/Measurement mode change
Measurement button + SET button	Memory mode: Calling up of stored values/setting of basic functions

The following symbols are used on the device, in the instructions for use, on the packaging, and on the type plate of the device:

⚠ WARNING

Indicates a potentially imminent danger. If it is not avoided, death or serious injury may result.

⚠ CAUTION

Indicates a potentially imminent danger. If it is not avoided, minor or minor injuries may result.

NOTE

Indicates a potentially harmful situation. If it is not avoided, the device or something in its vicinity may be damaged.



Product information

Note on important information



Follow the instructions

Read the instructions before starting work and/or operating equipment or machines



(Electrical) devices may not be disposed of with household waste



Manufacturer



Date of manufacture



CE marking

This product meets the requirements of the applicable European and national directives.



Labelling for identification of the packaging material

A = material abbreviation, B = material number:
1-7 = plastics, 20-22 = paper and cardboard



Device protected against solid foreign objects ≥ 12.5 mm and against dripping water at an angle



Direct current



Unique Device Identifier (UDI)

Marking for unique device identification



Item No.



Serial No.



Medical device



Model No.



Type BF applied part



Temperature range



Humidity range



Atmospheric pressure limits



Importer



Keep away from heat and sunlight



Protect from rain and keep dry



Authorised representative in the European Union

Available formats for the visually impaired:

The instructions for use available at www.aponorm.de in PDF format can be printed out in enlarged format.



www.aponorm.de

Intended use:

The thermometer is intended for intermittent measurement of human body temperature on the skin surface of the forehead. The thermometer is only intended for domestic/private use, and is not for commercial use.

Target group: The fever thermometer is suitable for newborns weighing 2.5 kg and over as well as infants, toddlers, and children under supervision, while self-use is possible for all individuals aged 12 and above. The thermometer may only be operated by persons understanding the instructions for use and who can operate the device safely. Certain groups of persons must be supervised by a person responsible for their safety while the thermometer is being used:

- Babies, toddlers, children
- Individuals with reduced capabilities (e.g. physical, mental, or sensory)

Indication:

Where an elevated body temperature is suspected. The user can use the thermometer to quickly and easily measure the human body temperature on the forehead. This allows the body temperature to be monitored.

Contraindications:

The measurement must not be performed during the healing phase following surgical procedures.

Clinical benefits:

The thermometer can be used to measure the body temperature on the forehead at home. Temperature monitoring is carried out in °C or °F with an accuracy of one decimal place.

Scope of delivery:

Check the scope of delivery for external damage of the cardboard packaging and completeness of the contents. Before use, ensure that the device and accessories and/or replacement parts have no visible damage and that all packaging material has been removed. If in doubt, do not use it and feel free to contact the [aponorm®](mailto:aponorm@customer-service.de) customer service via the product website www.aponorm.de/kontakt or the specialist retailer from whom you purchased the device.

- 1 x [aponorm®](http://www.aponorm.de) Contact Free 4 forehead thermometer
- 2 x 1.5 V AAA (LR03) batteries
- 1 x instructions for use + warranty card
- 1 x storage bag

Dear Customer,

Thank you for choosing the aponorm® Contact Free 4 forehead thermometer!

This forehead thermometer is a high-quality device for measuring body temperature for home use. It has been clinically tested and verified according to international standards.

This thermometer has been specially developed for safe use on the forehead. The device can be used to measure an individual's body temperature by detecting the intensity of infrared light emitted from the forehead. It converts the measured heat into a temperature value that can be read on the display. When used correctly, your temperature will be determined quickly and accurately.

Should you have any questions, problems or require spare parts,

please contact the aponorm® customer service department at any time via the product website www.aponorm.de/kontakt or the specialist retailer from whom you purchased the device. On the product website you will also find a lot of other useful information about your device.



www.aponorm.de/
kontakt

We wish you all the best for your health - your aponorm® team!

Currency of the instructions for use: 04/02/2025

Table of contents

1. Warnings and safety notes	10	6. Troubleshooting	25
2. Information on fever	15	7. Warranty	27
3. Operation	17	8. Disposal	29
- Inserting/changing batteries	17	9. Technical specifications	31
- Basic settings (unit, time & date) ..	18	10. Reporting of serious incidents	33
4. Application	20		
- Body mode	20		
- Object mode	21		
- Display measurements	22		
5. Cleaning & care	23		
- Clean device	23		
- Store device	24		

1. Warnings and safety instructions

General notes

⚠ WARNING

- Keep children away from the packaging material. There is a risk of suffocation.
- Keep the device away from children, pets, and vermin.
- Children must not use the device. Medical devices are not toys.
- The measurements taken by you are for your information only - they do not replace a medical examination! Discuss your measurements with your doctor and under no circumstances base your own medical decisions on them (e.g. regarding the dosage of medication)!
- If you notice any unusual temperature results, please contact your doctor immediately.
- The manufacturer assumes no liability for damage caused by improper or incorrect use.



- Store the device in the room where the measurement will be taken for at least 30 minutes beforehand.
- Wait approx. 30 minutes after exercising, taking a bath, or being outdoors before taking a measurement on the forehead.
- The thermometer must only be used to measure human body temperature.



- The thermometer is only intended for the measuring points on the human body specified in the instructions for use.
- The thermometer is calibrated during manufacture. If the thermometer is used according to the instructions for use, no further calibration is required.
- Check that the measuring sensor is intact before each use. If it is damaged, please contact the specialist retailer where you purchased the device.
- If you should have any questions about the use of the device, please contact customer service via the product website www.aponorm.de/kontakt at aponorm® or the specialist retailer from whom you purchased the device.

Risks for the device

NOTE

- Repairs may only be carried out by authorised service centres. Failure to respect this will void the warranty.
- Protect the device from impacts, moisture, dirt, strong temperature fluctuations, and direct sunlight.

- Do not expose the device to liquids. The device is not waterproof. Avoid any direct contact with water or other liquids.
- Portable and mobile HF communications equipment can affect the device.
- Do not touch the measuring sensor with your finger.
- Clean the thermometer after each use.

Notes on handling batteries

⚠ WARNING

- **Risk of explosion! Fire hazard!** Failure to observe the following points may result in personal injury or overheating, leakage, venting, breakage, explosion, or fire with the battery.
- This device contains non-rechargeable batteries. These must not be charged.
- Do not throw batteries into fire.
- Never charge, force discharge, heat, disassemble, open, crush, deform, encapsulate, or modify batteries.

- Never short-circuit batteries or battery compartment contacts.
- Protect the batteries from direct sunlight, rain, heat, and water.
- Exposure of batteries to an environment with extremely high temperatures or extremely low air pressure may result in explosion or leakage of flammable liquids and gases.
- Dispose of defective and depleted batteries immediately and properly (see section Disposal).
- Do not use modified or damaged batteries.
- Always select the correct battery type.
- Always insert batteries correctly and pay attention to the correct polarity (+/-).
- Never mix batteries of different manufacturers, capacities (new and used), sizes, and types within one device.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- If your skin or eyes come into contact with fluid from a battery cell, wash the affected areas with water and seek medical assistance.

- **Swallowing hazard!** Keep batteries out of reach of children. If batteries are swallowed, seek medical attention immediately. Swallowing can lead to burns, severe internal injuries and death.
- Never allow children to replace batteries without adult supervision.

⚠ CAUTION

- Store batteries in well-ventilated, dry, and cool rooms in a non-conductive container in which the batteries cannot short-circuit each other or be short-circuited by other metal objects.
- Keep batteries clean and dry.
- Keep batteries away from water.
- Remove the batteries from the battery compartment if the device is not to be used for an extended period of time.

NOTE

- Do not use rechargeable batteries.

Notes on electromagnetic compatibility

⚠ CAUTION

- The device is suitable for use in all environments listed in these instructions for use, including the domestic environment.
- In the presence of electromagnetic interference, the device may only be usable to a limited extent. As a result, error messages or a failure of the display/device may occur, for example.
- Use of this device directly adjacent to other devices, or stacked with other devices should be avoided because it could result in faulty operation. If, however, it becomes necessary to use the device as described above, both this device and any other devices should be monitored to ensure they are functioning properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to increased electromagnetic emissions or reduced electromagnetic immunity of the device and result in improper operation.
- Keep portable RF communication devices (including peripherals such as antenna cables or external antennas) at least 30 cm away from all parts

of the device, including all cables supplied.

- Non-compliance may lead to a reduction in the performance characteristics of the device.

2. Information about fever

There are many places where fever can be measured. The measured body temperature always depends on the measurement method (oral, axillary, rectal, on the forehead, in the ear). In general, do not compare measurement results from different measurement methods with each other, as the normal temperature range varies depending on the measurement location.

For this reason, contrary to widespread opinion, there is no universally recognised normal temperature. The normal temperature range is also influenced by other factors:

- Individual, person-dependent metabolism
- Age: The "normal" body temperature is higher in infants and small children than it is in adults. In children, higher temperature fluctua-

tions occur more quickly and more frequently and can also subside more quickly. As we get older, our "normal" body temperature drops.

- Clothing (babies: nappies as well for rectal measurement)
- Ambient temperature
- Time of day: Body temperature is lower in the morning and increases over the course of the day towards the evening.
- Activities: Physical activity and, with less influence, mental activity (e.g. stress) increase body temperature.

IMPORTANT: To determine your **personal normal temperature range**, ideally measure your body temperature several times on different days at the same time of day, **when you are feeling well**. In the event of illness, you can use these values as a basis for your fever rating.

NOTE

- **Rectal measurement** is considered the **most reliable method** for detecting fever, as it is very unlikely to be affected by user error.

If you are unsure about other measurement methods, you can reassure yourself by taking a rectal measurement. For **infants**, doctors recommend measuring **only rectally during the first 6 months**, as all other measurement methods can lead to inconsistent results.

You can also find more helpful information in our **7 golden rules for measuring fever**:





www.aponorm.de/
fiebermessregeln

3. Operation

Insert the batteries

1. Remove the battery compartment cover ④.
2. Insert the two batteries into the battery compartment with the poles facing the correct direction.
3. Only use batteries of the following type: 2 x 1.5 V AAA (LR03) batteries. Do not use rechargeable batteries.
4. Then close the battery compartment by fitting the cover and letting it click into place.


Replacing the batteries

When the batteries are low, this symbol will appear . Temperature measurements are still possible, but the batteries must be replaced soon. When  appears in the display, the batteries must be replaced. If the batteries are too low, the thermometer switches off automatically.

To replace the batteries, follow the previous steps under "Inserting the batteries" and also observe the following steps::

- Remove the depleted batteries from the battery compartment.
- Insert the new batteries. Make sure that the batteries are inserted with poles correctly aligned.

Switching on the thermometer

Turn the measuring sensor ① to the appropriate position to switch on the thermometer. After a short self-test and a short beep, the thermometer is ready to measure the temperature on the forehead. The device always starts in body mode .

Setting the basic function

Set the basic settings when using the thermometer for the first time. You can make the following basic settings: Temperature unit (celsius/

fahrenheit), date and time. With the thermometer switched off, press and hold the SET button for ③ 3 seconds to enter the setting mode.

Temperature unit

1. Press the measure button ②, to select the desired **unit**.
2. Confirm the unit by pressing the SET button ③. The time format display appears.

Setting the time format

After confirming the temperature unit, the time format appears automatically.

1. The device can display the time in either **12-hour** format (AM/PM) or in **24-hour**-format (24:00). Press the Measure button ② to select the format.
2. Confirm the format by pressing the SET button ③. The hour display starts to flash automatically.

Setting the time

After confirming the time format, the hour display starts to flash automatically.

1. Press the Measure button ②, to move forward one **hour**. Press the measure button ② repeatedly until the correct hour is displayed.
2. Confirm the hour setting by pressing the SET button ③. The minute display starts to flash automatically.

Follow the instructions under 1. & 2. to set the **minutes**.

Setting the date

After confirming the minutes, the **year** starts to flash automatically.

1. Press the measure button ② to move forward one year. Press the measure button ② repeatedly until the correct year is displayed.
2. Confirm the year setting by pressing the SET button ③. The month display starts to flash automatically.

Follow the instructions under 1. & 2. to set the **month** and the **day**.

Sound activation/deactivation

You can activate/deactivate sound on the device. If the sound is activated, a beep sounds after each measurement.

- To activate the sound, briefly press the SET button ③ when the device is switched on. A loudspeaker symbol appears on the display.
- To deactivate the sound, press the SET button ③ again. A changed loudspeaker symbol appears on the display.

4. Application

Body mode

Please consider that

- Physical activity, generally increased perspiration on the forehead, the use of vasoconstrictive medications, and skin irritations can distort the measurement result.
- the forehead must be free of sweat, hair, cosmetics, and dirt. If the forehead area is covered, clean it and then wait 10 minutes to improve the accuracy of the measurement.

- Keep the thermometer still during the measurement. Moving the thermometer during the measurement may cause inaccurate measurements.

NOTE

- We recommend that you take three measurements on the forehead. If the three measurements differ, take the highest measurement.

1. Turn the measurement sensor ① to the appropriate position to switch on the thermometer. The display is activated and shows all the segments. After the self-test, the adjacent figure appears on the display together with a beep. The measurement can now be started.

2. Point the thermometer at the centre of the forehead at a distance of less than 5 cm and then press the measurement button ②. If the distance is too great, "FAR" will appear on the display. The measurement starts automatically as soon as the measurement distance is correct. The measured value appears on the display and the end of the measurement is signalled with a short beep (if sound is activated).



3. Read the temperature on the display. In addition to the temperature, the colour display also indicates the fever level.


Colour	Temperature range	Fever indicator
Green	34.0 - 37.4 °C (93.2 - 99.3 °F)	Normal/no fever
Orange	37.5 - 38.5 °C (99.5 - 101.3 °F)	Increased temperature
Red	38.6 - 43.0 °C (101.5 - 109.4 °F)	High fever

If the temperature is within the fever range (red display), an acoustic fever indicator also sounds (4 x beeps).

4. The device switches off automatically after 60 seconds.
5. After each use, clean the device as outlined in the "Cleaning and Care" section.

Object mode

1. Turn the measurement sensor ① to the appropriate position to switch on the thermometer.

- Press the SET button ③ for 3 seconds to switch between Body and Object mode. After switching to Measurement mode , you can measure the object temperature.
- Orientate the thermometer at a distance of less than 5 cm from the centre of the object to be measured or the liquid (never immerse it in liquid!) and then press the measure button ②.
- The measured value appears together with a green backlight in the display and the end of the measurement is signalled with a short beep (if sound is activated).
- The device switches off automatically after 60 seconds.
- Clean the device after each use as described in the "Cleaning and care" section.

Displaying measurements

The device automatically stores the last 30 measurement results. If the 30 memory spaces are exceeded, the oldest value is deleted. The memory can be called up as follows:

- With the device switched on, briefly press the measure button ② and the SET button ③ simultaneously. The most recent measurement is dis-

played. The date and time, the memory space number, the symbol for the measurement mode, and the temperature are all displayed alternately.

- Each additional press of the SET button ③ displays the next stored measurement. The most recent measurement stored in the memory is marked with the number 1 and the oldest with the number 30.
- The user can ② press the measure button to take new measurements. The device switches off automatically after 60 seconds.

5. Cleaning and care

Cleaning of the device

NOTE

- The measuring sensor is the most sensitive part of the thermometer. Handle the measuring sensor with great care when cleaning.
- Do not use aggressive cleaning agents. Antiseptics, warm water, thinners, and benzene are not permitted.
- Always follow all the safety instructions for the user and the device ("Warnings & Safety Instructions").

Clean the measuring sensor after each use. Wipe the measuring sensor and the thermometer with a soft cloth that has been slightly moistened with 75% isopropyl alcohol solution. Do not use abrasive cleaners. Allow the device to dry for at least 10 minutes after cleaning before measuring the temperature.

NOTE

- Use only isopropyl alcohol to clean the measuring sensor. Clean the display and the outside of the thermometer with a soft, dry cloth.
- Under no circumstances should water be allowed to enter the device.
- Do not use the device again until it is completely dry.

Storing the device

NOTE



- The device must not be stored or used at a too high or too low temperature or humidity ("Technical Data"), in sunlight, together with electricity, or in dusty places. Otherwise measurement inaccuracies may occur.



- For planned long-term storage, remove the batteries.



6. Troubleshooting

If the forehead temperature reading is implausible or unrealistic, check that the forehead is free of sweat and cosmetics and that the forehead is not covered by hair. Follow the instructions for use and measure again.

Display	Cause	Remedy
Er 1	The thermometer is not working properly.	Remove the batteries, wait one minute, and then switch the device on again. If the message appears again, contact aponorm® Customer Service or the retailer from whom you purchased the device for assistance.

Display	Cause	Remedy
	The ambient temperature is outside the range of 5.0 °C to 40.0 °C (41.0 °F to 104.0 °F).	Leave the thermometer in a room with a temperature between 5.0 °C and 40.0 °C (41.0 °F to 104.0 °F) for at least 30 minutes.
	The measured temperature is higher than 1) Body mode: 43.0 °C (109.4 °F) 2) Object mode: 100.0 °C (212.0 °F)	Only operate the thermometer within the specified temperature ranges. Clean the measuring sensor if necessary. If the error message appears repeatedly, contact aponorm® customer service or the specialist retailer from whom you purchased the device.

Display	Cause	Remedy
	The measured temperature is lower than 1) Body mode: 34.0 °C (93.2 °F) 2) Object mode: 0 °C (32 °F)	Only operate the thermometer within the specified temperature ranges. Clean the measuring sensor if necessary. If the error message appears repeatedly, contact aponorm® customer service or the specialist retailer from whom you purchased the device.
	The thermometer is working properly.	Use the thermometer as usual.

Display	Cause	Remedy
	If the battery indicator flashes, this means that the batteries are low, but you can continue with the measurement.	The thermometer is measuring correctly, but the batteries need to be replaced soon.
	The thermometer does not work because the batteries are low.	Insert two new AAA alkaline batteries.

7. Warranty

We, WEPA Apothekenbedarf GmbH & Co. KG, voluntarily grant end customers who are consumers a 5-year warranty for all aponorm® Contact Free 4 forehead thermometers purchased in Germany and Austria in accordance with the following provisions:

Object and scope of the warranty

The warranty applies to all aponorm® Contact Free 4 forehead thermometers purchased by customers as consumers in Germany and Austria. WEPA Apothekenbedarf GmbH & Co. KG guarantees customers that the product is free from material and workmanship defects. If, despite proper handling of the product, a defect occurs within the warranty period of 5 years, WEPA Apothekenbedarf GmbH & Co. KG will repair the product free of charge (excluding transport costs for sending the device) or replace it completely at its own discretion.

Exclusions from the guarantee

Claims under this warranty are not valid in the event of damage caused by the customer or third parties, as a result e.g. of a fall, accident or improper handling. The same applies to damage caused by leaking batteries or

failure to follow the instructions for use. A warranty coverage by WEPA Apothekenbedarf GmbH & Co. KG is also excluded if a defect/damage is attributable to improper repair or other interventions by third parties. The warranty does not cover wear parts, batteries, or accessories as specified in the instructions for use.

Term of the warranty

The warranty is valid for a period of 5 years from the date of purchase (warranty period). The date of the purchase receipt or the warranty card filled out by the dealer with the date of purchase is decisive. The warranty period is not extended by any provision of services under this warranty, especially as it relates to repair or replacement of the device. The warranty period does not start over again in these cases either.

Assertion of warranty claims

The customer may exercise their warranty rights by presenting the defective product and proof of purchase, or a warranty card completed by the dealer, either to WEPA Apothekenbedarf GmbH & Co. KG or to the dealer from whom the product was purchased, within the warranty period. If a warranty claim is valid, WEPA Apothekenbedarf GmbH & Co. KG will either repair the product free of charge (excluding transport costs

for sending in the device) or replace it entirely at its discretion. If the check reveals that there is no valid warranty claim, the product can be repaired by WEPA Apothekenbedarf GmbH & Co. KG at the customer's expense. In this case, WEPA Apothekenbedarf GmbH & Co. KG shall inform the customer of the expected costs by means of a cost estimate before carrying out the repair. The customer may then choose to have the product repaired at their own expense or request the return of the unrepaired product.

Further rights of the customer

WEPA Apothekenbedarf GmbH & Co. KG provides the warranty under the stated conditions. In addition, the customer may be entitled to further legal rights. The customer's contractual and statutory rights, especially statutory warranty rights against the product's seller, remain unaffected by this warranty. The warranty is subject to the laws of the Federal Republic of Germany.

8. Disposal

In the interest of protecting the environment, do not dispose of the device with household waste at the end of its service life.

Disposal can be carried out at appropriate collection points in your country.



Dispose of the device in accordance with EC Directive - WEEE (Waste Electrical and Electronic Equipment).

■ If you have any questions, please contact the local authorities responsible for waste disposal.

- The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points, or at electronics retailers. You are legally obliged to dispose of the batteries.
- You will find these symbols on batteries containing harmful substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury.



Information for consumers on the disposal and scrapping of old appliances in Germany

Owners of old appliances can hand them in free of charge within the framework of the options set up and provided by public disposal authorities for the return or collection of old appliances, to ensure their proper disposal. The return is required by law.

In addition, returns can be made to distributors under certain conditions. Pursuant to the German Electrical and Electronic Equipment Act (ElektroG), the following distributors are obliged to take back old equipment free of charge:

- Specialised electrical stores with a sales area for electrical and electronic equipment of at least 400 square metres.
- Grocery retailers with a total sales area of at least 800 square metres that offer electrical and electronic devices several times a year or on a permanent basis and place them on the market.
- Online retail businesses, where the obligation for one-for-one take-back in private households applies only to heat exchangers (refrigerators/freezers, air conditioners, etc.), display devices, and large appliances. For the one-for-one take-back of lamps, small appliances, and small IT and telecommunications equipment, as well as zero-for-one take-back, online retailers must provide return options within a reasonable distance from the end consumer.

These distributors are obliged,

- When handing over a new electrical or electronic device to an end user, an old device of the same type from the end user, which essentially serves the same function as the new device, must be taken back free of charge either at the place of handover or in its immediate vicinity, and
- to take back, free of charge, and at the end-user's request, old appliances that do not exceed 25 centimetres in their external dimensions in the retail store or in the immediate vicinity thereof; the take-back may not be linked to the purchase of an electrical or electronic appliance and is limited to three old appliances per type of appliance. You can find out about local return options from your retailer.

If the products contain batteries and rechargeable batteries or lamps that can be removed from the old device without destruction, these must be removed before disposal and disposed of separately as batteries or lamps.

If the old electrical or electronic equipment contains personal data, you are responsible for deleting it before returning it.

9. Technical data

If the device is used outside of its specifications, proper functioning cannot be guaranteed!

The accuracy of this thermometer has been carefully tested and developed for a long service life.

When using the device in medical applications, the applicable national regulations determine whether metrological checks must be carried out using appropriate means.

We reserve the right to make technical changes to improve and further develop the product.

Basic functions	Forehead temperature measurement
Measurement units	Celsius (°C) or Fahrenheit (°F)
Mode of functioning	Adjustment mode
Conditions of usage	5.0 °C to 40.0 °C (41.0 °F - 104.0 °F) at 15-85% relative humidity and 700-1060 hPa ambient pressure

Storage and transportation conditions	-25.0 °C to 55.0 °C (-13.0 °F - 131.0 °F) at a relative humidity of 15-95%
Measuring distance	1-5 cm distance to measurement point
Minimum measurement duration	1 second
Measuring range:	Body mode: 34.0 °C to 43.0 °C (93.2 °F - 109.4 °F) Object mode: 0.0 °C to 100.0 °C (32.0 °F - 212.0 °F)
Accuracy of measurement	Body mode: ± 0.2 °C (0.4 °F) from 35.5 °C to 42.0 °C (95.9 °F - 107.6 °F) outside this measuring range ± 0.3 °C (± 0.5 °F) Object mode: ± 4% or ± 2.0 °C (± 4.0 °F)
Clinical measurement accuracy (Body mode)	Clinical systematic error: -0.23 °C (-0.41 °F) Clinical reproducibility: 0.1 °C (0.18 °F) Conformity limit: 0.67 °C (1.21 °F)
Body mode memory function	Automatically stores the last 30 measurements.

Display	LCD-display
Display resolution	0.1 °C/°F
Power-saving features	The device switches off automatically after 60 seconds.
Dimensions (W x D x H)	approx. 153.9 x 41.4 x 25.9 mm
Weight	approx. 102 g (with batteries)
Batteries	2 x 1.5V AAA (LR03) operating time for approx. 6,000 measurements
IP protection	IP22
Expected service life	5 years

The serial number can be found on the device or in the battery compartment.

This device complies with the European standard EN 60601-1-2 (Group 1, Class B, in compliance with CISPR11, IEC61000-4-2, IEC61000-4-3, IEC61000-4-8) and is subject to particular precautions with regard to electromagnetic compatibility.

Please note that portable and mobile RF communication equipment can interfere with the functioning of this device.

10. Reporting serious incidents

If serious incidents occur during the use of the product, such as a deterioration in health or unusual incidents, please report this to WEPA Apothekenbedarf GmbH & Co. KG. You may also report the incident to the competent authority of the respective EU Member State. In the Federal Republic of Germany this is:

Federal Institute for Drugs and Medical Devices
Pharmacovigilance Department
Kurt-Georg-Kiesinger Avenue 3
D-53175 Bonn
Germany
www.bfarm.de

Incidents are considered serious if they have directly or indirectly caused, could cause, or might cause harm

- a) the death of a person,
- b) a temporary or permanent serious deterioration in a person's state of health, or
- c) a serious risk to public health has occurred, could have occurred, or could occur.

If you have any other questions about the product or need technical support, please use the contact form on the product website www.aponorm.de/kontakt to contact aponorm® customer service.



Notes

Product features



Distance checking

Only measures when the correct distance to the centre of the forehead is detected.



Self-measurements

Facilitates self-measurement thanks to the rotatable measuring sensor.



Measurement memory

Stores the last 30 measurements.



Non-contact measurement

Allows measurement while the child is sleeping without disturbing them.



3 in 1 measuring range

Body, surface and ambient temperature measurable.



Accuracy of measurement

± 0.2 °C in the range 35.5 - 42.0 °C



Colour display

Illuminated display with fever indicator enables measurements in the dark.



Signal tone can be deactivated

Silent mode enables interference-free measurement.

**CLINICALLY +
PROVEN**

5 year
warranty



DIE APOTHEKENMARKE

Distributed by:
WEPA Apothekenbedarf
GmbH & Co KG, D-56204
Hillscheid, Germany
www.wepa-apothekenbedarf.de



JOYTECH Healthcare Co., Ltd.
No. 365, Wuzhou Road, 311100
Hangzhou, Zhejiang Province, PCR



Shanghai
International Holding Corp.
GmbH (Europe), Eiffestraße 80,
20537 Hamburg, Germany

CE 0123



Koviva GmbH
Bahnstraße 16,
40212 Düsseldorf,
Germany

Ordering No.: 046712



PZN -13659829