
Ohr-Thermometer comfort

technology by
microlife

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
Ear Thermometer comfort

technology by
microlife

EN The apornorm[®] by microlife Ear Thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this thermometer can provide a stable, heat-interference-free reading with each measurement. The instrument performs a self-test every time it is switched on to always guarantee the specified accuracy of measurements. The apornorm[®] by microlife Ear Thermometer is intended for the periodic measurement and monitoring of human body temperature in the home. It is intended for use on people of all ages. **This thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.** Please read through these instructions carefully so that you understand all functions and safety information.

The Advantages of this Ear Thermometer

Multiple Uses (Wide Range Measurement)

This thermometer offers a wide measurement range feature from 0 °C to 100.0 °C (32.0 °F to 212.0 °F), meaning the unit can be used as an ear thermometer to measure body temperature or it can be used to measure surface temperature of the following:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Probe Cover Free

This thermometer is more user-friendly and more cost-effective since a probe cover is not required.

Measurement in 1 second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

Accurate and reliable

The unique probe assembly construction, incorporating an advanced infrared sensor, ensures that each measurement is accurate and reliable.

Gentle and Easy to Use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child, thereby causing no disruption.
- This thermometer is quick, therefore child-friendly.

Auto-Display Memory

The last reading is automatically displayed for 2 seconds when the unit is switched ON.

Multiple Reading Recall

Users will be able to recall the last 12 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and Hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe can be done with an alcohol-moistened cotton tissue, making this thermometer completely hygienic for use by the whole family.

Fever Alarm

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature above 37.5 °C.

Important Safety Instructions

- This instrument may only be used for the purposes described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- **Never immerse this instrument in water or other liquids. For cleaning please follow the instructions in the «Cleaning and Disinfecting» section.**
- Do not use the instrument if you think it is damaged or notice anything unusual.
- Never open the instrument.
- Earwax in the ear canal may cause a lower temperature reading. Therefore it is important to ensure the subject's ear canal is clean.
- This instrument comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Protect it from:
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
- If the instrument is not going to be used for a prolonged period the battery should be removed.
- Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

WARNING: Use of this instrument is not intended as a substitute for consultation with your physician. This instrument is NOT waterproof so do not immerse in liquids.

How this Ear Thermometer measures Temperature

This thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) ensures the most accurate ear temperature. Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of fever.

To avoid an inaccurate measurement:

1. Switch on the thermometer by pressing the «ON/OFF» button.
2. After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.
3. Place the probe firmly into the ear canal, press the «START» button and keep the probe in the ear until the thermometer beeps to signal the completion of the measurement.

Control Displays and Symbols

- **All segments displayed:** Press the «ON/OFF» button to turn on the unit; all segments will be shown for 2 seconds.
- **Memory:** The last reading will automatically be shown on the display for 2 seconds.
- **Ready for measurement:** When the unit is ready for measurement, the «°C» or «°F» icon will keep flashing.
- **Measurement complete:** The reading will be shown on the display with the «°C» or «°F» icon; the unit is ready for the next measurement, when the «°C» or «°F» icon is flashing again.
- **Out-of-ear temperature indication:** A crossed-ear icon will appear on the display if the reading falls outside the range 32.0 ~ 42.2 °C (89.6 ~ 108.0 °F).
- **Low battery indication:** When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery.

Directions for Use

1. Press the «ON/OFF» button. The display is activated to show all segments for 2 seconds.
2. The last measurement reading will automatically be shown on the display for 2 seconds with the «M» icon.
3. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measurement.
4. Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.
 - For children under 1 year; pull the ear straight back.
 - For children of 1 year to adult; pull the ear up and back.

5. While gently pulling the ear, insert the probe snugly into the ear canal (max. 1 sec.) and immediately press the «START» button. Release the button and wait for the beep sound. This is the indication that confirms the end of measurement.
6. Remove the thermometer from the ear canal. The display shows the measured temperature.

NOTE:

- **To ensure accurate readings, wait at least 30 sec. after 3-5 continuous measurements.**
- Accumulation of ear wax on the probe can result in less accurate temperature readings or cross infection between users. **Therefore, it is essential that the probe is clean before each measurement.** For cleaning, follow the instructions in the «Cleaning and Disinfecting» section.
- **After cleaning the measuring sensor with alcohol, wait 5 minutes before taking the next measurement,** in order to allow the thermometer to reach its operating reference temperature.
- 10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature above 37.5 °C.
- For an infant, it is best to have the child lying flat with his/her head sideways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
- Always take the temperature in the same ear, since the temperature readings may be different from ear to ear.
- Wait for a few minutes to take the ear temperature after sleeping.
- In the following situations it is recommended that three temperatures in the same ear be taken and the highest one taken as the reading:
 1. New born infants in the first 100 days.
 2. Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 3. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.
 4. If the measurement is surprisingly low.

Changing between Celsius and Fahrenheit

This thermometer can display temperature measurements in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn «OFF» the unit, press and hold the «START» button for 5 seconds; after 5 seconds, the current measurement scale («°C» or «°F» icon) will be flashing on the display. Change the measurement scale between °C and °F by pressing the «START» button. When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the «ready for measuring» mode.

How to recall 12 readings in Memory Mode

This thermometer can recall the last 12 readings.

- **Recall mode:** Press the «START» button to enter Recall mode when power is off. The memory icon «M» will flash.
- **Reading 1 – the last reading:** Press and release the «START» button to recall the last reading. Display alone with memory icon.
- **Reading 12 – readings in succession:** Press and release the «START» button consecutively to recall the last 12 readings in succession.

Cleaning and Disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the probe lens and the display.

Battery Replacement

This instrument is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the display. Remove the battery cover by sliding it in the direction shown. Remove the battery and replace with a new one.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

| | Technical Specifications |
|-------------------------------|--|
| Type | Ear Thermometer comfort |
| Measurement range | 0 °C to 100.0 °C (32.0 °F to 212.0 °F) |
| Resolution | 0.1 °C / °F |
| Measurement accuracy | Laboratory: ±0.2 °C, 32.0 ~ 42.2 °C (±0.4 °F, 89.6 ~ 108.0 °F) ±1 °C, 0 ~ 31.9 °C, 42.3 ~ 100.0 °C (±2 °F, 32.0 ~ 89.5 °F, 108.1 ~ 212.0 °F) |
| Display | Liquid Crystal Display, 4 digits plus special icons |
| Acoustic | <ul style="list-style-type: none"> • The unit is turned ON and ready for the measurement: 1 short beep • Complete the measurement: 1 long beep • System error or malfunction: 3 short beeps • Fever alarm: 10 short beeps |
| Memory | <ul style="list-style-type: none"> • Auto-Display the last measured temperature • 12 readings recall in the Memory Mode |
| Backlight | <ul style="list-style-type: none"> • The display light will be GREEN for 4 seconds, when the unit is turned ON. • The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C (99.4 °F). • The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C (99.4 °F). |
| Operating temperature | 5 °C to 40 °C (41.0 °F to 104 °F) |
| Storage temperature | -25 °C to +55 °C (-13 °F to 131 °F) |
| Automatic Switch-off | Approx. 1 minute after last measurement has been taken. |
| Battery | CR2032 Battery (X1) 3V – at least 1000 measurements |
| Dimensions | 140 x 47 x 15 mm |
| Weight: | 59 g (with battery), 56 g (w/o battery) |
| Reference to standards | Complies with ASTM E1985 requirements EN 12470-5, ASTM E 1985, IEC 60601-1, IEC 60601-1-2 |

www.apornorm.de

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.apornorm.de.

Technical alterations reserved.

According to the Medical Product User Act a biennial technical inspection is recommended for professional users. Please observe the applicable disposal regulations.

Guarantee

We grant you a **5 year guarantee** after the date of purchase. Any damage caused by improper handling shall not be covered by the warranty. Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

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WEPA Apothekenbedarf GmbH & Co. KG