
Thermometer flexible Kupfer

technology by
microlife

DE

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EN

Introduction of Antimicrobial Copper Digital Thermometer

This Digital Antimicrobial Medical Thermometer provides quick and highly accurate readings over the human body temperature range. At the same time, it reduces the microbial flora, by the dispersion of contagious microorganisms, thus providing to be safer for the end user. Please read the instructions thoroughly before use. The quality of the thermometer has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive) Annex I essential requirements and applied harmonized standards. EN 12470-3: 2000 / A1: 2009 Clinical thermometers- Part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.

Antimicrobial Copper Properties

Surfaces made or covered by special copper alloys, has strong antimicrobial properties against a wide variety of microorganisms*. Copper alloys emit copper ions Cu⁺ that whilst in contact with microbes and bacteria rupture their cellular membranes, thus destroying these microorganisms. This activity reduces the microbial flora on the coated area and – due to the "halo phenomenon" – simultaneously causes a drastic reduction in pathogens on the remaining body of the thermometer. Thermometers containing parts with antimicrobial copper alloys drastically reduce microbial flora, dispersing contagious microorganisms, thus proving to be safer for the end user.

* Efstathiou A. Panos «The Role of Antimicrobial Copper Surfaces in Reducing Healthcare associated Infections», European Infectious Disease, Volume 5, Issue 2, Autumn 2011

Important safety instructions

- The instrument may be used only for measuring body temperature!
- The minimum measurement time until the beep is heard must be maintained without exception!
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m. In cases where you suspect this to be unavoidable, please verify if the device is working properly before use.
- Protect the instrument from impact and dropping!
- Avoid bending the tip more than 45°!
- Avoid ambient temperatures above 60 °C. NEVER boil the instrument!
- The manufacturer recommends to verify the accuracy by an authorised laboratory every 2 years.

⚠ Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Antimicrobial Thermometer Use

Recommendation: Before use, keep the antimicrobial thermometer out of any physical contact, for at least 2 hours, this utilizes the antimicrobial copper properties, as described in the above section "Antimicrobial Copper Properties". The use of Antimicrobial copper is a supplement to and not a substitute for standard infectious control practices; Users must continue to follow all current infection control and cleaning practices.

Turning on the Thermometer

To turn on the thermometer, press the ON/OFF button ①; a short beep signals «thermometer ON». A display test is performed. All segments should be displayed. The last measurement reading will be shown on the display ② automatically for 2 seconds with the «M» icon. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «C» appear at the display field ②. The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

Using the Thermometer

Choose the preferred measuring method. When taking a measurement, the current temperature is continuously displayed and the «C» symbol flashes. If the beep is heard 10 times and the «C» is no longer flashing, the predictive end-temperature has been determined and the thermometer can be read now. 10 short beeps will sound when the temperature is higher than 37.5 °C in order to alert the patient that he/she may have fever. Reference: Oral temperature. To achieve comparable results allow a 1 minute interval time between measurements. To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button ①. Otherwise the thermometer will automatically turn off after about 10 minutes.

Measuring methods / Normal body temperature

In the mouth (oral) / 35.5 - 37.5 °C

Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor ④ must be in good contact with the tissue. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. **Approx. measuring time: 10 seconds!**

In the anus (rectal) / 36.6 - 38.0 °C

This is the most reliable measuring method, and is especially suitable for infants and small children. Carefully insert the measuring sensor ④ of the thermometer 2 to 3 cm into the anal aperture. **Approx. measuring time: 10 seconds!**

In the armpit (axillary) / 34.7 - 37.3 °C

To receive more reliable results aponorm® by microlife recommends measuring temperature orally or rectally.

Cleaning and Disinfecting

Clean the thermometer with a soft, dry cloth or with a cotton tissue moistened with Isopropyl alcohol (70%). Don't let the thermometer come into contact with any chemical thinner. **Please never immerse into liquids!**

Battery Replacement

When the «▼» symbol (upside-down triangle) appears at the display field, the battery is flat and needs replacing. To replace the battery remove the battery compartment cover ③ from the thermometer. Insert the new battery with the + at the top. Make sure you have a battery of the same type to hand. Batteries can be purchased at any electrical store.

Type	Technical specifications
Measurement range	Predictive maximum thermometer 32,0 °C to 42,9 °C Temp. < 32,0 °C: display «L» for low (too low) Temp. > 42,9 °C: display «H» for high (too high)
Measurement accuracy	± 0,1 °C between 34 °C and 42 °C
Operating conditions	10 - 40 °C; 15-95 % relative maximum humidity
Storage conditions	-25 - +60 °C; 15-95 % relative maximum humidity
Battery	1,5/1,55 V; SR41
Battery lifetime	approx. 2700 measurements (using a new battery)
IP Class	IP22
Storage conditions	temperature: -25 °C to +60 °C, humidity: 15 % to 95 % (noncondensing)
Reference to standards	EN 12470-3, clinical thermometers; ASTM E1112; IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11
Expected service life	5 years or 10000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

Guarantee

We grant you a warranty of 5 years after the date of purchase. Any damage caused by improper handling shall not be covered by the warranty. Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

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WEPA Apothekenbedarf GmbH & Co. KG, 56204 Hillscheid (Germany),